

APLW 004

Mobile unit

for optimum loading of large-volume laboratory glassware, up to 50 l.

EAN: 4002516115427 / Material number: 11104510 / Old Material Number: 69300400D



- Equipped with 1 nozzle

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6111	•
Product category	
Injector basket	•
Module for large-volume laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Flasks [number]	1
Dimensions and weight	
External dimensions, net height in mm	588
External dimensions, net width in mm	535
External dimensions, net depth in mm	582
External dimensions, gross height in mm	760
External dimensions, gross width in mm	625
External dimensions, gross depth in mm	715
Net weight in kg	4,68
Gross weight in kg	9,83