


APLW 008
 Mobile unit
 for optimum loading of 3 pipette cartridges.



- For lower level
- Dimensions of pipette cartridge: W 84 x H 384 x D 84 mm

EAN: 4002516115465 / Material number: 11104550 / Old Material Number: 69300800D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6111	•
Product category	
Mobile unit for pipettes	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Area of application	
Reprocessing of pipettes/phials	•
Dimensions and weight	
External dimensions, net height in mm	456
External dimensions, net width in mm	179
External dimensions, net depth in mm	366
External dimensions, gross height in mm	355
External dimensions, gross width in mm	565
External dimensions, gross depth in mm	610
Net weight in kg	7,64
Gross weight in kg	9,64