

APLW 015
Upper basket
with 16 unoccupied injector slots.



- Only suitable for use on upper levels

EAN: 4002516115571 / Material number: 11104650 / Old Material Number: 69301500D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Upper basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	16
Flasks [number]	16
Erlenmeyer flasks [number]	16
Round flasks [number]	16
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	532
External dimensions, net depth in mm	560
External dimensions, gross height in mm	340
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	600
Net weight in kg	4,00
Gross weight in kg	6,60