

APLW 021

Lower basket

with 110 unoccupied injector slots for nozzles with Ø 2.5 mm.



- Only suitable for use on the lowest level

EAN: 4002516115762 / Material number: 11104740 / Old Material Number: 69302100D

| | |
|--|-------------------------------|
| Product classification | |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PLW 6011 | • |
| PLW 6111 | • |
| Product category | |
| Injector basket | • |
| Lower basket | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black, White |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Test tubes [number] | 110 |
| Phials [number] | 110 |
| Narrow-necked glassware [number] | 110 |
| Flasks [number] | 110 |
| Dimensions and weight | |
| External dimensions, net height in mm | 575 |
| External dimensions, net width in mm | 504 |
| External dimensions, net depth in mm | 450 |
| External dimensions, gross height in mm | 355 |
| External dimensions, gross width in mm | 562 |
| External dimensions, gross depth in mm | 615 |
| Net weight in kg | 5,91 |
| Gross weight in kg | 8,35 |