

APLW 024/1  
Lower basket  
equipped with 42 APLW 052 injector nozzles.



- Only suitable for use on the lowest level

EAN: 4002516115809 / Material number: 11104780 / Old Material Number: 69302401D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in mm	160
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	352
External dimensions, gross width in mm	562
External dimensions, gross depth in mm	615
Net weight in kg	5,22
Gross weight in kg	8,83
Companion pack	
APLW 052 injector nozzle, Ø 4 mm, length 175 mm [number]	42