

APLW 024/3

Lower basket

equipped with 30 APLW 066 injector nozzles and 12 APLW 067 injector nozzles.



**Image
in progress**

- Only suitable for use on the lowest level

EAN: 4002516115823 / Material number: 11104800 / Old Material Number: 69302403D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector basket	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel, Black, White
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	42
Flasks [number]	42
Erlenmeyer flasks [number]	42
Round flasks [number]	42
Dimensions and weight	
External dimensions, net height in mm	220
External dimensions, net width in mm	535
External dimensions, net depth in mm	575
External dimensions, gross height in mm	352
External dimensions, gross width in mm	562
External dimensions, gross depth in mm	615
Net weight in kg	5,22
Gross weight in kg	13,87
Companion pack	
APLW 066 injector nozzle, Ø 6 mm, length 175 mm [number]	30
APLW 067 injector nozzle, Ø 6 mm, length 225 mm [number]	12