


APLW 033
 Top frame
 to create an even surface.



- Suitable for the APLW 000
- Reduction in maximum load height by 50 mm

EAN: 4002516115953 / Material number: 11104930 / Old Material Number: 69303300D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other inserts	•
Retaining frame	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	444
External dimensions, net width in mm	473
External dimensions, net depth in mm	60
External dimensions, gross height in mm	190
External dimensions, gross width in mm	435
External dimensions, gross depth in mm	512
Net weight in kg	1,22
Gross weight in kg	2,85