

APLW 036

Lid

Height 45 mm, for securing delicate wash items.



EAN: 4002516115984 / Material number: 11104960 / Old Material Number: 69303600D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Other inserts	•
Cover	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	248
External dimensions, net width in mm	248
External dimensions, net depth in mm	45
External dimensions, gross height in mm	278
External dimensions, gross width in mm	57
External dimensions, gross depth in mm	285
Net weight in kg	0,39
Gross weight in kg	0,53