

APLW 045

Injector nozzle

with plastic support, conical, Ø 2.5, length 80 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 4 mm

EAN: 4002516116073 / Material number: 11105050 / Old Material Number: 69304500D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	15
External dimensions, net width in mm	15
External dimensions, net depth in mm	90
External dimensions, gross height in mm	40
External dimensions, gross width in mm	65
External dimensions, gross depth in mm	118
Net weight in kg	0,01
Gross weight in kg	0,02