

APLW 048

Injector nozzle

with plastic support, star-shaped, Ø 2.5, length 80 mm.

EAN: 4002516116097 / Material number: 11105080 / Old Material Number: 69304800D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 4 mm

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	32
External dimensions, net width in mm	32
External dimensions, net depth in mm	90
External dimensions, gross height in mm	5
External dimensions, gross width in mm	110
External dimensions, gross depth in mm	167
Net weight in kg	0,00
Gross weight in kg	0,00