

APLW 053  
 Injector nozzle  
 with plastic support, star-shaped, Ø 6, length 115 mm.

EAN: 4002516116134 / Material number: 11105130 / Old Material Number: 69305300D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 10 mm

|  |                          |
|--|--------------------------|
| <b>Product classification</b>                  |                          |
| Large-capacity washer-disinfectors, laboratory | •                        |
| <b>Machine affiliation</b>                     |                          |
| PLW 6011                                       | •                        |
| PLW 6111                                       | •                        |
| <b>Product category</b>                        |                          |
| Injector nozzle                                | •                        |
| <b>Product properties</b>                      |                          |
| Material                                       | Stainless steel, Plastic |
| Colour   | Blue, Stainless Steel    |
| <b>Area of application</b>                     |                          |
| Reprocessing of laboratory glassware           | •                        |
| <b>Dimensions and weight</b>                   |                          |
| External dimensions, net height in mm          | 76                       |
| External dimensions, net width in mm           | 76                       |
| External dimensions, net depth in mm           | 123                      |
| External dimensions, gross height in mm        | 55                       |
| External dimensions, gross width in mm         | 180                      |
| External dimensions, gross depth in mm         | 250                      |
| Net weight in kg                               | 0,00                     |
| Gross weight in kg                             | 0,01                     |