

APLW 056

Injector nozzle

with plastic support, star-shaped with clamp, Ø 6, length 225 mm.

EAN: 4002516116219 / Material number: 11105160 / Old Material Number: 69305600D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 10 mm

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	75
External dimensions, net width in mm	75
External dimensions, net depth in mm	232
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Net weight in kg	0,02
Gross weight in kg	0,03