

APLW 064

Injector nozzle

with plastic support and bottle neck holder, star-shaped, Ø 6, length 225 mm.

EAN: 4002516116301 / Material number: 11105240 / Old Material Number: 69306400D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 10 mm
- With locking clamp

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	76
External dimensions, net width in mm	76
External dimensions, net depth in mm	231
External dimensions, gross height in mm	55
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	300
Net weight in kg	0,04
Gross weight in kg	0,10