

APLW 065

Injector nozzle

with plastic support and bottle neck holder, star-shaped, Ø 6, length 275 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 10 mm
- With locking clamp

EAN: 4002516116318 / Material number: 11105250 / Old Material Number: 69306500D

<b>Product classification</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 6011	•
PLW 6111	•
<b>Product category</b>	
Injector nozzle	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Blue, Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	75
External dimensions, net width in mm	75
External dimensions, net depth in mm	287
External dimensions, gross height in mm	55
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	300
Net weight in kg	0,04
Gross weight in kg	0,11