

APLW 068  
 Injector nozzle  
 with holder, Ø 6, length 275 mm.



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Including nozzle cap, Ø 10 mm

EAN: 4002516116349 / Material number: 11105280 / Old Material Number: 69306800D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product properties	
Material	Wire, rilsan-coated, Stainless steel
Colour	Blue, Stainless Steel, White
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	70
External dimensions, net width in mm	85
External dimensions, net depth in mm	287
External dimensions, gross height in mm	55
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	300
Net weight in kg	0,09
Gross weight in kg	0,15