

APLW 071

Injector nozzle

w. holder for meas. cylinders up to Ø 105 mm & injector nozzles Ø 8, l. 320 mm.

EAN: 4002516116370 / Material number: 11105310 / Old Material Number: 69307100D



- Suitable for standard mobile injector units and injector baskets
- Suitable e.g. for measuring cylinders, flasks and test tubes
- Contact surface at tip of nozzle, Ø 17 mm

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Injector nozzle	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Grey
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	120
External dimensions, net width in mm	120
External dimensions, net depth in mm	330
External dimensions, gross height in mm	125
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	350
Net weight in kg	0,24
Gross weight in kg	0,38