

APLW 086

Bottle neck holder

Ø 45 mm, for injector nozzles Ø 6 mm.



Image
in progress

EAN: 4002516116646 / Material number: 11105460 / Old Material Number: 69308600D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Other accessories	•
Product properties	
Material	Plastic
Colour	White
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	33
External dimensions, net width in mm	33
External dimensions, net depth in mm	70
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0,05
Gross weight in kg	0,05