

APLW 088

Bracket

for measuring cylinders up to Ø 105 mm and injector nozzles Ø 8 mm.



EAN: 4002516116660 / Material number: 11105480 / Old Material Number: 69308800D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Holder	•
Other accessories	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	12
External dimensions, net width in mm	120
External dimensions, net depth in mm	300
External dimensions, gross height in mm	108
External dimensions, gross width in mm	117
External dimensions, gross depth in mm	305
Net weight in kg	0,14
Gross weight in kg	0,23