

APLW 090

Holder

for pipettes with maximum 14 mm and minimum 4 mm diameter.



EAN: 4002516116684 / Material number: 11105500 / Old Material Number: 69309000D

Product classification	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 6011	•
PLW 6111	•
Product category	
Holder	•
Other accessories	•
Product properties	
Material	Plastic
Colour	Grey, Black
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	20
External dimensions, net width in mm	20
External dimensions, net depth in mm	20
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0,00
Gross weight in kg	0,00