


A 302/3

Module

for optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Equipped with 16 A 844 nozzles and 16 A 845 nozzles
- Part of the EasyLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth

EAN: 4002516188988 / Material number: 11116960 / Old Material Number: 69630203D

| | |
|---|--------------------------|
| Product category | |
| Washer-disinfectors, laboratory | • |
| Large-capacity washer-disinfectors, laboratory | • |
| Machine affiliation | |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| Product category | |
| Insert/module for narrow-necked glassware | • |
| Insert for various flasks | • |
| Product properties | |
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |
| Part of the EasyLoad system | • |
| Area of application | |
| Reprocessing of laboratory glassware | • |
| Capacity | |
| Narrow-necked glassware [number] | 32 |
| Dimensions and weight | |
| External dimensions, net height in mm | 143 |
| External dimensions, net width in mm | 238 |
| External dimensions, net depth in mm | 479 |
| External dimensions, gross height in mm | 135 |
| External dimensions, gross width in mm | 290 |
| External dimensions, gross depth in mm | 590 |
| Net weight in kg | 1,49 |
| Gross weight in kg | 2,67 |
| Companion pack | |
| A 844, injector nozzle Ø 2.5 mm, length 80 mm [number] | 16 |
| A 845, injector nozzle Ø 2.5 mm, length 125 mm [number] | 16 |