

AK 12/1
Insert
for loading different utensils.



- Suitable for beakers, funnels, breast milk pumps, for example
- With flexible handles

EAN: 4002516343714 / Material number: 11131690 / Old Material Number: 69501204D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8592	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert/module for beakers	•
Insert/module for wide-necked glassware	•
Other inserts	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	122
External dimensions, net width in mm	224
External dimensions, net depth in mm	434
External dimensions, gross height in mm	125
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	450
Net weight in kg	0,65
Gross weight in kg	1,07