

A 857

Conversion kit

For using lab carts with the G 7827/28 and PG 8527/28 machines in the PG 86xx.



- Use of existing Miele load carriers in the PG 86xx series
- For all G 7827, G 7828, PG 8527 and PG 8528 load carriers

EAN: 4002516201946 / Material number: 11235960 / Old Material Number: 69685700D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Conversion kit	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	60
External dimensions, gross width in mm	110
External dimensions, gross depth in mm	405
Net weight in kg	0,33
Gross weight in kg	0,43