


A 11/2

Insert

For the gentle cleaning of delicate load items by interrupting the water jet.



- With perforations 7 x 7 mm
- For use in upper and lower basket

EAN: 4002516340119 / Material number: 11239130 / Old Material Number: 69501103D

Product classification	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8505	•
Product category	
Underlay	•
Retaining frame	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, net height in mm	1
External dimensions, net width in mm	476
External dimensions, net depth in mm	443
External dimensions, gross height in mm	70
External dimensions, gross width in mm	505
External dimensions, gross depth in mm	430
Net weight in kg	0,71
Gross weight in kg	1,08