

LDTT 86

Collection bottle

With frame for collecting drained liquids from the TT 86.



**Image  
in progress**

- Installed underneath the drip tray

EAN: 4002516271864 / Material number: 11293740 / Old Material Number: 69200120D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Storage/transport system	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Dimensions and weight</b>	
External dimensions, net height in mm	121
External dimensions, net width in mm	130
External dimensions, net depth in mm	147
External dimensions, gross height in mm	140
External dimensions, gross width in mm	190
External dimensions, gross depth in mm	295
Net weight in kg	0,41
Gross weight in kg	0,56