

A 320

Insert

For the gentle cleaning of delicate load items by interrupting the water jet.



- With perforations 7 x 7 mm

EAN: 4002516237082 / Material number: 11295510 / Old Material Number: 69632000D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Retaining frame	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	10
External dimensions, net width in mm	237
External dimensions, net depth in mm	536
External dimensions, gross height in mm	60
External dimensions, gross width in mm	245
External dimensions, gross depth in mm	560
Net weight in kg	0,47
Gross weight in kg	0,63