

DOS K 85/1 flex

Dispenser module for 5 litre canister

with level indicator for the automatic dispensing of liquid media.

EAN: 4002516365303 / Material number: 11498930 / Old Material Number: 69747053USA



- Dispenser monitoring – high process security as per EN ISO 15883
- For dispensing liquid cleaning agent/neutralising agent
- 214 mm siphon for 5 l media canisters with flexible cap
- Siphon conversion kit for 10–30 l media canisters available

<b>Product classification</b>	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8583	•
PG 8591	•
PG 8592	•
PG 8593	•
<b>Product category</b>	
Dispenser system	•
<b>Product properties</b>	
Production rate in ml/60 s	55
Length of suction lance in mm	214
<b>Standard electrical connection</b>	
Number of phases	1
Neutral conductor	N
AC/DC	AC
Voltage in V	200-240
Frequency in Hz	60
Total rated load in kW	0,01
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	160
External dimensions, gross width in mm	198
External dimensions, gross depth in mm	478
Net weight in kg	2,22
Gross weight in kg	2,44
<b>Standards, certs, symbols</b>	
CE	•