

RF MF27/28

Adapter

for using the MF 27/28 transfer trolley with the PLW 8615/16/17 model series.



- Retrofit release element for the MF 27/28 transfer trolley
- Installed by a service technician
- Two are required for the double-door PLW 86165 machines

EAN: 4002516352471 / Material number: 11499280 / Old Material Number: 69249201D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Adapter	•
<b>Product properties</b>	
Material	Plastic
Colour	Black
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	12
External dimensions, gross width in mm	120
External dimensions, gross depth in mm	170
Net weight in kg	0,25
Gross weight in kg	0,25