

DKF 30-R

Reactivatable Longlife AirClean charcoal filter
for Miele recirculation cooker hoods DAD 4840 and DAD 4940.



EAN: 4002516382836 / Material number: 11662830 / Old Material Number: 28996418EU1

- Air so fresh and odour-free that you can feel it in your kitchen
- Active charcoal for effective neutralisation of unpleasant odours
- Can easily be regenerated in the oven every six months
- Service life 3 years
- Contains: 2 filters without frames
- Special construction type for low flow resistance