

Boost CX1 Cat & Dog PowerLine  
 Bagless cylinder vacuum cleaners  
 for the highest hygiene demands, in a compact design.

EAN: 4002516384649 / Material number: 11666830 /  
 Old Material Number: 41NCF031GB

<b>General information</b>	
Model identifier	Boost CX1 Cat & Dog PowerLine
Manufacturer material number	11666830
EAN number	4002516384649
<b>Construction type</b>	
Bagless vacuum cleaner	•
Product brand	Boost CX1
Motor type	PowerLine
Product type	SNCF0
<b>Appliance type and design</b>	
Appliance colour	Obsidian black
<b>Cleaning</b>	
Vortex Technology	•
<b>User convenience</b>	
Electronic suction power control	4-stage rotary selector
Ergonomic carrying handle	•
Ergonomic handle	•
Suction tube	Comfort telescopic tube
Park system for storing	on one side
<b>Mobility</b>	
TrackDrive	•
Operating radius in m	10
Cable length in m	6,5
<b>Filter system</b>	
Exhaust filter	HEPA AirClean filter
Dustbag volume in l	1,0
<b>Technical data</b>	
Power	890
Volume in dB	78,0
Width in mm	400
Height in mm	280
Depth in mm	280
Weight in kg	7,8
<b>Accessories included</b>	
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush	•
Floorhead SBD 365-3	•
Turbobrush STB 305-3 (DOUBLE; LOCKED)	•
Cleaning tool	No

- Powerful cleaning performance thanks to vortex technology
- Always stays on track thanks to TrackDrive
- Maximum air hygiene with HEPA filter: \*Ideal for allergy sufferers
- with additional TurboBrush – ideal for pet lovers
- Extremely convenient height adjustment thanks to a Comfort telescopic tube.

\* According to EN 60312-1