

A 865
Holder
for holding multi-functional and 3F syringes perfectly.



- For use in A 105/1 upper basket and A 315 module
- Compatible with syringes from several well-known manufacturers

EAN: 4002516413752 / Material number: 11709860 / Old Material Number: 69686500D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8581	•
PG 8591	•
Product category	
Holder	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	121
External dimensions, net width in mm	28
External dimensions, gross height in mm	108
External dimensions, gross width in mm	117
External dimensions, gross depth in mm	305
Net weight in kg	0,30
Gross weight in kg	0,93