

A 830/1

Insert

ideal holder for separable mouth specula.



- For use in insert E 337/1 and E 802/1
- 3 compart. Ø 8 mm, 14 compart. 10 x 19 mm, 2 compart. 8 x 19 mm
- Mouth specula handles can be placed in compartments of inserts
- Suitable for mouth specula ≥ 3

EAN: 4002516429203 / Material number: 11729480 / Old Material Number: 68683001D

Product category	
Thermal disinfectors	•
Machine affiliation	
PG 8528	•
G 7831	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Insert/module for dental instruments	•
Product properties	
Material	Plastic
Colour	Black
Number of compartments	16
Area of application	
Reprocessing of dental instruments	•
Dimensions and weight	
External dimensions, net height in mm	14
External dimensions, net width in mm	49
External dimensions, net depth in mm	159
External dimensions, gross height in mm	67
External dimensions, gross width in mm	170
External dimensions, gross depth in mm	230
Net weight in kg	0,02
Gross weight in kg	0,02