


APLW 215

Injector module

for the optimum loading of large-volume plastic bottles up to 50 l.

EAN: 4002516525981 / Material number: 11853940 / Old Material Number: 69621500D



- For use on level 1
- Equipped with a special injector nozzle
- Equipped with a holding frame
- Max. bottle diameter: 380 mm
- Bottle height: 560 to 620 mm
- Bottle neck diameter: min. 30 mm internal, max. 80 mm external

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
Product category	
Other inserts/other modules	•
Module for large-volume laboratory glassware	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Large-volume laboratory glassware [number]	1
Dimensions and weight	
External dimensions, net height in mm	565
External dimensions, net width in mm	469
External dimensions, net depth in mm	526
External dimensions, gross height in mm	590
External dimensions, gross width in mm	520
External dimensions, gross depth in mm	530
Net weight in kg	3,84
Gross weight in kg	5,14