


APLW 152
Basket
for various inserts.



- For use on level 1
- Vertical clearance dependent on other baskets used
- Clearances: width 507 mm, depth 545 mm

EAN: 4002516526001 / Material number: 11854310 / Old Material Number: 69615200D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 7111	•
Product category	
Carrier/basket carrier	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Dimensions and weight	
External dimensions, net height in mm	180
External dimensions, net width in mm	559
External dimensions, net depth in mm	568
External dimensions, gross height in mm	180
External dimensions, gross width in mm	570
External dimensions, gross depth in mm	625
Net weight in kg	2,33
Gross weight in kg	3,41