

B 4847 FashionMaster  
Steam ironing system  
with display and steamer for perfect ironing results and convenience.



- Uniform steam distribution – Miele honeycomb soleplate
- Professional results thanks to 4.0 bar steam pressure
- Convenient and easy to set up/store – 1-2 lift system
- Top-class ironing results thanks to active ironing board
- Convenient smoothing of hanging textiles – Steamer

EAN: 4002516527220 / Material number: 11978160 / Old Material Number: 13484721AS1

<b>Design</b>	
Iron colour	Anthracite/lotus white
Ironing board cover colour	Tayberry red
<b>Perfect ironing results</b>	
Honeycomb soleplate	•
Steam pressure in bar	4,0
Steam quantity in g/min	100
Vertical steam	•
Steam starting from temperature level	2
Steam burst	•
Continuous steam	•
Steamer	•
Active ironing table (inflation and suction function)	•
Suction function	•
2-stage fan	•
Non-stick soleplate	•
Proven virus hygiene	•
<b>Time-saving ironing</b>	
Water container volume in l	1,25
Removable water container	•
Permanently refillable water container	•
Multi-layer ironing	•
<b>User convenience</b>	
All-in-one	•
1-2 Lift system	•
Variable height adjustment with gas pressure stay, in mm	830-1020
Castors gentle on the floor	•
Ironing board with Comfort Zone	•
Fan safety switch-on/switch-off	•
Fan memory function	•
Choice of language	•
Guided first-time use	•
Guided descaling	•
Automatic descaling	•
Automatic rinsing	•
Storage compartment for iron	•
Storage compartment for mains cable and accessories	•
<b>Safety</b>	
CoolDown function	•
AutoOff function	•
<b>Technical data</b>	
Dimensions in mm (width)	470
Dimensions in mm (height)	970
Dimensions in mm (depth)	1480
Dimensions in mm (width), folded	470
Dimensions in mm (height), folded	1280
Dimensions in mm (depth), folded	370
Frequency in Hz	50-60
<b>Accessories supplied</b>	
Non-stick soleplate	•

B 4847 FashionMaster  
 Steam ironing system  
 with display and steamer for perfect ironing results and convenience.



EAN: 4002516527220 / Material number: 11978160 / Old Material Number: 13484721AS1

Steamer	•
Premium ironing table cover	•
Steam hose holder	•
Test strips for determining water hardness	•
Descaling tablets	•

**Available display languages**

Available display languages via MultiLingua

bahasa malaysia, dansk, deutsch, english (AU), english (GB), english (US), español, français, hrvatski, italiano, magyar, nederlands, norsk, polski, portugês, русский, română, slovenčina, slovenščina, srpski, suomi, svenska, türkçe, українська, čeština, ελληνικά,