

PKF-RO/1

Active charcoal filter

To protect the membranes of RO systems.

EAN: 4002516691655 / Material number: 12295400 / Old Material Number: 67809802D



Product classification	
Tank dishwashers	•
Machine affiliation	
PTD 701	•
PTD 702	•
Product category	
Water treatment	•
Dimensions and weight	
External dimensions, gross height in mm	95
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	285
Net weight in kg	0,49
Gross weight in kg	0,58

- Removes substances affecting taste and smell, cloudiness, and Cl
- Flow rate range: 0–300 l/h
- Working pressure: 1–6 bar. Max. drop in pressure: 40 kPa
- Water temperature: 4–35 °C
- Maximum lifetime: 1 year or after filtering 40 m³ of water