


APTD 221  
Basket set  
for mixed load items in bistros and cafés.



- Base basket incl. 3 inserts for plates, glasses and cutlery
- Highly versatile
- Base basket dimensions: 500 x 500 mm

EAN: 4002516717812 / Material number: 12305290 / Old Material Number: 67022101D

<b>Product category</b>	
Tank dishwashers	•
Throughfeed tank dishwashers	•
<b>Machine affiliation</b>	
PTD 901	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•
<b>Product category</b>	
Basket for mixed tableware	•
<b>Product properties</b>	
Material	Plastic
Colour	Grey
Number of rows	3
<b>Dimensions and weight</b>	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	225
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	1,92
Gross weight in kg	2,96