

A 100

Upper basket

For the optimum loading of injector modules or inserts.



- Upper basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Loading height with module with plastic supports 190 mm
- Clearances: width 489, length 440 mm

EAN: 4002516841296 / Material number: 12642640 / Old Material Number: 69610002D

Product category

Washer-disinfectors, laboratory

•

Cleaning machines for industrial parts cleaning

•

Machine affiliation

PLW 8636 [LAB / LAB MON]

•

PG 8504

•

PG 8583

•

PG 8583 CD

•

PG 8593

•

PLW 8604

•

PLW 8683

•

PLW 8683 CD

•

PLW 8693

•

Product category

Basket with module connection

•

Upper basket

•

Product properties

Material

Stainless steel, Plastic

Colour

Stainless Steel, Black

Area of application

Reprocessing of pipettes/phials

•

Reprocessing of laboratory glassware

•

Dimensions and weight

External dimensions, net height in mm

141

External dimensions, net width in mm

525

External dimensions, net depth in mm

528

External dimensions, gross height in mm

170

External dimensions, gross width in mm

555

External dimensions, gross depth in mm

565

Net weight in kg

2,92

Gross weight in kg

3,56