

Guard L1 Allergy
Cylinder vacuum cleaner
powerful | HEPA AirClean filter | ideal for allergy sufferers



- 2-in-1 upholstery and crevice nozzle for easy cleaning without switching over
- Ideal power for every cleaning task: 4 selectable modes
- Premium design concept with stylish applications
- Individually selectable height and ergonomic Comfort handle
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal for allergy sufferers thanks to HEPA AirClean filter
- Durable: tested for the equivalent of 20 years* use
- ComfortFit: developed for the highest hygiene standards

* During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More:

miele.com/20years

EAN: 4002516832706 / Material number: 12652230 /
Old Material Number: 41UFF033GB

Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard L1
Motor type	PowerLine
Product type	SUFFO
Construction type and design	
Appliance colour	Brilliant white
User convenience	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Suction tube	Suction tube X(HES)
Park system for pauses during vacuuming	•
Comfort cable rewind	•
Comfort handle 2.0	•
Mobility	
Castors with steel shaft, gentle on floors	•
4-point manoeuvrability	•
Operating radius in m	12
Cable length in m	8,5
Filter system	
Vacuum cleaner bag type	HyClean Pure TU
Dustbag volume in l	3,5
Dustbag change indicator	•
Exhaust filter	HEPA AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
Technical data	
Power	890
Volume in dB	74,0
Width in mm	500
Height in mm	239
Depth in mm	28,80
Weight in kg	6,6
Accessories included	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush	•
Floorhead SBD 365-3	•
HyClean Pure TU vacuum cleaner bag	•
Voucher for HyClean Pure TU	•