

Guard M1 Flex

Cylinder vacuum cleaner

powerful suction | integrated accessories | including flexible crevice nozzle



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years** use
- Height-adjustable EasySlide telescopic tube
- Ideal for hard-to-reach places: flexible crevice nozzle
- Mode button: simply adjust the suction power to your needs
- Practical compartment with three-piece accessory set for added convenience
- Stylish design concept and illuminated LED icons
- Range extra, fewer socket changes: 11 m operating radius

* During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More:

miele.com/20years

EAN: 4002516832980 / Material number: 12652550 /
Old Material Number: 41VWF034GB

General information	
Model identifier	Guard M1 Flex
Manufacturer material number	12652550
EAN number	4002516832980
Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard M1
Motor type	PowerLine
Product type	SVWFO
Construction type and design	
Colour	Nordic blue
User convenience	
Electronic suction power control	MODE button
Ergonomic carrying handle	•
Ergonomic handle	•
Suction tube	HES mid-size telescopic tube
Park system for pauses during vacuuming	•
Park system for storing	on both sides
Comfort cable rewind	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	11
Cable length in m	7,5
Filter system	
Vacuum cleaner bag type	HyClean Pure CO
Dustbag volume in l	4,5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
Furniture protection	•
Technical data	
Power	890
Volume in dB	77,0
Width in mm	496
Height in mm	227
Depth in mm	284
Weight in kg	6,7
Accessories included	
Accessory storage	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Extra-long, flexible crevice nozzle SFD 20	•
Dusting brush	•

Guard M1 Flex

Cylinder vacuum cleaner

powerful suction | integrated accessories | including flexible crevice nozzle



EAN: 4002516832980 / Material number: 12652550 /

Old Material Number: 41VWF034GB

Floorhead SBD 365-3

•

HyClean Pure CO vacuum cleaner bag

•

Voucher for HyClean Pure CO

•