

Guard S1 Junior

Cylinder vacuum cleaner

high suction power | light and compact | in modern trend colours | easy to use



- ComfortFit: developed for the highest hygiene standards
- Ergonomic handle: for comfortable handling
- 99.99 %* dust retention capacity with AirClean filter
- Durable: tested for the equivalent of 20 years*** use
- 10 m operating radius: effortless cleaning of larger areas
- Easy handling thanks to lightweight and compact design

* In accordance with IEC 62885-2:2021

** During the development phase of the advertised Guard vacuum cleaner series, Miele tested models and core components at various power levels for 750 hours of use (= 37.5 hours p.a.). More:

miele.com/20years

EAN: 4002516833727 / Material number: 12652860 /

Old Material Number: 41TAF030GB

Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard S1
Motor type	PowerLine
Product type	STAFO
Construction type and design	
Appliance colour	Terra red
User convenience	
Electronic suction power control	4-stage rotary selector
Suction tube	EasyLock sectioned suction tube
Park system for pauses during vacuuming	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	10
Cable length in m	7,5
Filter system	
Vacuum cleaner bag type	HyClean Pure TU
Dustbag volume in l	3,5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
Furniture protection	•
Technical data	
Module depth in mm	24,50
Power	890
Volume in dB	76,0
Width in mm	425
Height in mm	230
Weight in kg	6,0
Accessories included	
Accessory storage	On VarioClip
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush	•
Floorhead SBD 355-3	•
HyClean Pure TU vacuum cleaner bag	•
Voucher for HyClean Pure TU	•