

A 860
Holder rack
for A 300/2 and A 300/3 injector modules.



- Arrangement: 2 x 4 fields
- For use on A 300/2 and A 300/3 injector modules
- For the optimum positioning of laboratory glassware
- Part of the EasyLoad system

EAN: 4002516855118 / Material number: 12665410 / Old Material Number: 69686001D

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 7111	•
Product category	
Other inserts/other modules	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	8
Compartment width in mm	114
Compartment length in mm	114
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Dimensions and weight	
External dimensions, gross height in mm	170
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	490
Net weight in kg	0,48
Gross weight in kg	0,88