

E 139/1

Insert

For the optimum loading of test tubes up to 12 x 200 mm.



- Suitable for test tubes up to 12 x 200 mm
- Divided into 6 surfaces
- Mesh width 9 x 9 mm
- For use in lower basket
- Including lid A 13

EAN: 4002516854036 / Material number: 12667390 / Old Material Number: 69513903D

Product category

Washer-disinfectors, laboratory

•

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

G 7825

•

PG 8527

•

PG 8535

•

PG 8536

•

PLW 8636 [LAB / LAB MON]

•

PG 8504

•

PG 8583

•

PG 8583 CD

•

PG 8593

•

PLW 8615

•

PLW 8616

•

PLW 8617

•

PLW 8683

•

PLW 8683 CD

•

PLW 8693

•

PLW 6011

•

PLW 6111

•

PLW 7111

•

PLW 8505

•

Product category

Insert/module for test tubes

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Area of application

Reprocessing of test tubes

•

Reprocessing of laboratory glassware

•

Capacity

Test tubes [number]

200

Dimensions and weight

External dimensions, gross height in mm

245

External dimensions, gross width in mm

245

External dimensions, gross depth in mm

265

Net weight in kg

1,00

Gross weight in kg

1,21

Companion pack

A 13 lid [number]

1