


A 302/2

Injector module

for optimum loading of laboratory glassware – can be equipped with 32 nozzles.



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 32 injector nozzles
- Part of the EasyLoad system
- Nozzle spacing 60 mm across the width, 57 mm across the depth
- Particularly suitable for A 844 and A 845 (nozzle diameter 2.5 mm)

EAN: 4002516854098 / Material number: 12669210 /
Old Material Number: 69630204D

Product category

| | |
|---|---|
| Washer-disinfectors, laboratory | • |
| SlimLine laboratory washers | • |
| Cleaning machines for industrial parts cleaning | • |
| SlimLine laboratory washers for industrial parts cleaning | • |

Machine affiliation

| | |
|--------------------------|---|
| PLW 8636 [LAB / LAB MON] | • |
| PG 8504 | • |
| PG 8583 | • |
| PG 8583 CD | • |
| PG 8593 | • |
| PLW 8604 | • |
| PLW 8683 | • |
| PLW 8683 CD | • |
| PLW 8693 | • |
| PLW 7111 | • |

Product category

| | |
|---|---|
| Insert/module for narrow-necked glassware | • |
| Insert for various flasks | • |

Product properties

| | |
|-----------------------------|--------------------------|
| Material | Stainless steel, Plastic |
| Colour | Stainless Steel, Black |
| Part of the EasyLoad system | • |

Area of application

| | |
|--------------------------------------|---|
| Reprocessing of laboratory glassware | • |
|--------------------------------------|---|

Capacity

| | |
|----------------------------------|----|
| Narrow-necked glassware [number] | 32 |
|----------------------------------|----|

Dimensions and weight

| | |
|---|------|
| External dimensions, net height in mm | 72 |
| External dimensions, net width in mm | 195 |
| External dimensions, net depth in mm | 479 |
| External dimensions, gross height in mm | 90 |
| External dimensions, gross width in mm | 205 |
| External dimensions, gross depth in mm | 500 |
| Net weight in kg | 1,11 |
| Gross weight in kg | 1,96 |