

TU Allergy XL HyClean Pure
HyClean Pure Allergy XL pack
8 vacuum cleaner bags and 1 SF-HA 50 HEPA AirClean filter at a discount price.

EAN: 4002516925651 / Material number: 12694920 /
Old Material Number: 41996762INT



- Best filtration performance of 99.99 %* thanks to Miele AirClean system
- Exceptional service life thanks to patented 3D airflow guides
- HyClose automatic self-sealing collar for even greater hygiene
- ComfortFit – 1 of 11 patented details of the Miele HyClean® vacuum cleaner bag
- Maximum air hygiene with HEPA filter: **Ideal for allergy sufferers
- Suitable for all Guard L1 and Guard S1 models

* According to IEC 62885-2 when using Miele appliances fitted with AirClean filter system

** According to EN 60312-1