

CO Allergy XL HyClean Pure
HyClean Pure Allergy XL pack
8 vacuum cleaner bags and 1 SF-HA 50 HEPA AirClean filter at a discount price.

EAN: 4002516889366 / Material number: 12694930 /
Old Material Number: 41996763INT



- Best filtration performance of 99.99 %* thanks to Miele AirClean system
- Exceptional service life thanks to patented 3D airflow guides
- HyClose automatic self-sealing collar for even greater hygiene
- ComfortFit** – 1 of 11 patented details of Miele's HyClean® vacuum cleaner bag
- Maximum air hygiene with HEPA filter: ***Ideal for allergy sufferers
- Suitable for all Guard M1 models

* According to IEC 62885-2 when using Miele appliances fitted with AirClean filter system

** EP 2 744 383 B1

*** According to EN 60312-1