

CO XXL HyClean Pure

16 HyClean Pure CO vacuum cleaner bags at a special price

Miele AirClean system delivers the best in filtration performance at 99.99 %.

EAN: 4002516897521 / Material number: 12695230 /
Old Material Number: 41996765INT



- 80 % recycled material* and ultimate performance – unrivalled on the market
- Exceptional service life thanks to patented 3D airflow guides
- HyClose automatic self-sealing collar for even greater hygiene
- ComfortFit – 1 of 11 patented details of the Miele HyClean® vacuum cleaner bag
- Contains: 16 vacuum cleaner bags, 4 exhaust filters
- Suitable for all Guard M1 models

* with a tolerance of +/-5 %