

E 130  
Insert  
for the optimum loading of 10 trays.



- equipped with 11 holders, height 170 mm, distance apart 35 mm
- For use in lower basket

EAN: 4002511129375 / Material number: 02648950 / Old Material Number: 69513001D

<b>Product classification</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert for trays	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	35
Rail height in mm	170
<b>Area of application</b>	
Reprocessing of trays	•
<b>Capacity</b>	
Trays [quantity]	10
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	180
External dimensions, net width in mm	180
External dimensions, net depth in mm	445
External dimensions, gross height in mm	210
External dimensions, gross width in mm	210
External dimensions, gross depth in mm	470
Net weight in kg	0.76
Gross weight in kg	0.90