

E 363

Insert

for the optimum loading of different instruments.



- Mesh width 1 mm
- Equipped with lid
- For use in upper and lower basket

EAN: 4002511116283 / Material number: 03149790 / Old Material Number: 69536301D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	55
External dimensions, net width in mm	150
External dimensions, net depth in mm	225
External dimensions, gross height in mm	80
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	250
Net weight in kg	0.62
Gross weight in kg	0.71