

E 379

Insert

For holding different instruments.



- 0.8 mm wire mesh, 1.7 mm mesh width, 5 mm frames
- equipped with 2 handles
- For use in upper and lower basket

EAN: 4002511116450 / Material number: 03203500 / Old Material Number: 69537901D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Mesh tray	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	110
External dimensions, net width in mm	180
External dimensions, net depth in mm	445
External dimensions, gross height in mm	135
External dimensions, gross width in mm	235
External dimensions, gross depth in mm	485
Net weight in kg	0.72
Gross weight in kg	1.25