

E 378  
Insert  
For holding different instruments.



- 0.8 mm wire mesh, 1.7 mm mesh width, 5 mm frames
- equipped with 2 handles
- For use in lower basket

EAN: 4002511116443 / Material number: 03203530 / Old Material Number: 69537801D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
SlimLine thermal disinfectors, dental	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 7121	•
PWD 7122	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Area of application	
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of surgical instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in mm	110
External dimensions, net width in mm	460
External dimensions, net depth in mm	460
External dimensions, gross height in mm	165
External dimensions, gross width in mm	535

E 378  
Insert  
For holding different instruments.

EAN: 4002511116443 / Material number: 03203530 / Old Material Number: 69537801D

External dimensions, gross depth in mm	560
Net weight in kg	1.79
Gross weight in kg	2.25

